


# Gastrostomitube med silikon ballong



Produktnavn	Volum i ballong Min /Anbefalt /max	Antall i eske	Varern
AMTGastrostomitube 12FR	2 ml /2,5 ml / 3 ml	1	7-1210
AMTGastrostomitube 14FR	3 ml/ 4 ml / 5 ml	1	7-1410
AMTGastrostomitube 16FR	4 ml/ 6 ml / 8 ml	1	7-1610
AMTGastrostomitube 18FR	6 ml/ 8 ml / 10 ml	1	7-1810
AMTGastrostomitube 20FR	7 ml/ 10 ml / 15 ml	1	7-2010
AMTGastrostomitube 24FR	7 ml/ 10 ml / 15 ml	1	7-2410