

5th Generation Product Data Sheet

7L

Product Code:

5110XXXX

Please note the last four digits denote the specific aperture type and lid colour.

Our Sharpsafe containers are specifically designed to enable the safe disposal of clinical sharps and minimise the risk of associated injuries.

Our 7L is designed for low to medium use areas and ideal for acute and community settings.

Available with different colour lids for easily identifiable waste segregation. Products come with a range of compatible accessories for increased stability, security and safety.

Sharpsafe's new 5th Generation containers launched in 2022 and were the world's first full range of sustainable single use sharps containers. Engineered using low carbon footprint recycled material showing our commitment to championing change.



SUPPORTING
A CIRCULAR
ECONOMY




REDUCING
YOUR CARBON
FOOTPRINT



REDUCING VIRGIN
PLASTIC ENTERING
THE SUPPLY CHAIN



MADE USING LOCALLY
SOURCED LOW
CARBON FOOTPRINT
RECYCLED MATERIALS





| Product Certification | |
|-----------------------|---|
| ISO 23907-1:2019 | ✓ |
| UN3291 | ✓ |
| BSI Kitemark | ✓ |

| Site Certification | |
|--|---|
| ISO 9001 - Quality Management | ✓ |
| ISO 14001 - Environmental Management | ✓ |
| ISO 45001 - Occupational Health & Safety | ✓ |

| Product Details | |
|---------------------------|-------------------------|
| Size: | 7L |
| Usable Capacity: | 5.78L |
| Usable Volume: | 85% |
| Weight: | 401g/410g |
| Fill Line Height: | 308mm |
| Autoclavable: | Yes |
| Autoclavable Conditions: | 134°C |
| Product Dimensions (mm): | 380 x 175 x 175 (HxWxD) |
| Aperture Dimensions (mm): | 110 x 60 (LxW) |
| GTIN: | 5056271506173 |
| Lot/Batch Number: | Not Applicable |
| Expiry Date: | Not Applicable |

| Accessory Compatibility | | | |
|-------------------------|---|---------------------|---|
| Wall Bracket | ✓ | Trolley Bracket | ✓ |
| Stand | ✓ | Sharpsafe Trolley** | ✓ |
| NPSD Tray | ✗ | Pedestal Trolley** | ✓ |

| Logistic Data | |
|---------------------------------|-------------------------|
| Carton Quantity: | 50 |
| Carton Weight (when empty): | 1.6kg |
| Carton Weight (when full): | 23.5kg |
| Carton Dimensions (mm): | 680 x 345 x 810 (HxWxD) |
| Carton Volume: | 0.19m ³ |
| Storage Temperature Conditions: | Ambient |
| Storage Lower Temp Limit: | None |
| Storage Upper Temp Limit: | None |
| Sterile/Non Sterile: | Non Sterile |

| Features and Benefits | | | |
|----------------------------|-----|-------------------------------|---|
| In-Mould Label | ✓ | Latex Free | ✓ |
| Temporary Closure | ✓ | Low Carbon Footprint Material | ✓ |
| Permanent Closure | ✓ | Leak Resistance | ✓ |
| Four Click Assembly | ✓ | Disconnection Features | ✓ |
| Fill Line | ✓ | Safety Shield | ✓ |
| Protected Access** | n/a | Semi Translucent Lid | ✓ |
| Single Use | ✓ | Carry Handle | ✓ |
| Different Lid Variations** | ✓ | Puncture Resistant | ✓ |

| Pallet Information | |
|----------------------------------|-------------|
| Euro Pallet Carton Quantity: | 9 |
| Euro Pallet Dimensions (mm): | 1200 x 800 |
| Standard Pallet Carton Quantity: | 12 |
| Standard Pallet Dimensions (mm): | 1200 x 1000 |

| Manufacturing Information | |
|---------------------------|-----------|
| Made In: | Wales, UK |
| By: | Vernacare |

*Weight without safety shield/Weight including safety shield
 **Different product codes may apply
 ***When used with a Safe Clip Trolley Bracket Adaptor

FOR MORE INFORMATION:

✉ sharpsafesales@vernagroup.com

🌐 vernacare.com 📞 +44 (0)1495 235800