

# Sage Oral Care Petite System q4



**REF 6714-X**

## Specifications

Product name	Sage Oral Care Petite System q4
Product code	6714-X
Description	Daily kit to provide oral hygiene including tools for each 4-hour nursing shift
Manufacturer accreditation	ISO 13485:2003
Storage information	Do not store above 25° C
Expiration date	24 months from date of manufacture
Length of Suction Toothbrush	14.1 cm
Length of each bristle	8 mm
Length of Suction Swab	18.4 cm
Length of Applicator Swab	15.9 cm
Amount of solution in each burst packet	7 mL
Amount of Mouth Moisturizer	2 g
Case quantity	15

## Features and benefits

### Suction Toothbrush

Ultra-soft nylon bristles	<ul style="list-style-type: none"> <li>Gently removes dental plaque, debris, and oral secretions, known to harbor respiratory pathogens</li> </ul>
User-friendly thumb port	<ul style="list-style-type: none"> <li>Provides easy suction control</li> </ul>
Three suction ports	<ul style="list-style-type: none"> <li>Avoids clogging and keeps open path for debris</li> </ul>
Connects directly to standard suction lines	<ul style="list-style-type: none"> <li>Helps promote compliance</li> </ul>
Available with suction handle	<ul style="list-style-type: none"> <li>Provides variable suction control</li> </ul>
Angled handle	<ul style="list-style-type: none"> <li>Provides cleansing of the back of the oral cavity</li> </ul>

### Suction Swab

Suction Swab	<ul style="list-style-type: none"> <li>Helps remove debris and secretions between brushings</li> </ul>
Soft foam heads	<ul style="list-style-type: none"> <li>Gentle to fragile tissue</li> </ul>

## Alcohol-free Mouthwash ingredients

Water (aqua)  
 Sorbitol  
 Flavor (aroma)  
 Potassium Sorbate  
 Polysorbate 80  
 Polysorbate 20  
 Cetylpyridinium chloride  
 Citric acid  
 Blue 1 (CI 42090)

## Mouth Moisturizer ingredients

Water (aqua)  
 Cocos nucifera (coconut) oil  
 Xylitol  
 Flavor (aroma)  
 Cellulose gum  
 Tocopheryl acetate  
 Mentha viridis (spearmint) leaf oil  
 Polysorbate 20  
 Polysorbate 80  
 Potassium sorbate  
 Carbomer  
 Cetylpyridinium chloride

# Sage Oral Care Petite System q4

## Packaging information

### Contents:

- (2) packages of (1) Untreated Petite Suction Toothbrush, 7 mL Alcohol-Free Mouthwash, 2 g Mouth Moisturizer and (1) Applicator Swab.
- (4) packages of (1) Untreated Petite Suction Swab, 7 mL Alcohol-Free Mouthwash, 2 g Mouth Moisturizer and (1) Applicator Swab.

### Indications of use:

**Suction Toothbrush:** Aids in the removal of dental plaque, debris, and oral secretions.

**Suction Swab:** Aids in the removal of oral secretions and debris.

### Instructions:

**Suction Toothbrush:** Instructions: Burst solution packet with thumbs. Peel package open. Attach Suction Toothbrush to suction line. Brush teeth and oral cavity for 1 minute. To suction, place thumb over port. Rinse to clear tubing after use. Discard. Dispense Mouth Moisturizer onto Applicator Swab and apply to lips and oral cavity as needed

**Suction Swab:** Instructions: Burst solution packet with thumbs. Peel package open. Attach Suction Swab to suction line. Clean mouth for 1 minute. To suction, place thumb over port. Rinse to clear tubing after use. Discard. Dispense Mouth Moisturizer onto Applicator Swab and apply to lips and oral cavity as needed.

### Caution:

Do not allow patient to bite down on the oral care tool. Use a bite block if patient has altered levels of consciousness or cannot comprehend commands. Use caution with children and unresponsive individuals. Failure to follow these safety precautions may damage the device and present a choking/aspiration hazard. May not function as intended/potential risk of cross-contamination if device is reused.

Not made with natural rubber latex.

Suction Toothbrush and Suction Swab: MDD 93/42/EEC

Applicator Swab: MDD 93/42EEC

Alcohol-Free Mouthwash and Mouth Moisturizer: EC Cosmetic Regulation 1223/2009

### This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, Perox-A-Mint. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This document is not intended for distribution outside EU and EFTA.



Sage Products LLC  
3909 Three Oaks Road  
Cary, IL 60013 USA  
t: 0808 189 0456

[sage-products.co.uk](http://sage-products.co.uk)

24677B-UK  
Copyright © 2018 Stryker