

GO VOM® - Vomit Bag


An effective solution to biological waste management and containment



50 bags
per display

Bag : 6" x 15" (unfolded)

Display : 5,5" x 7,25" x 8"



- Hygienic and practical, this GO VOM® bag with absorbent pad is ideal for emergency interventions (hospital, medical clinic, ambulance), and while travelling by: plane, car, boat etc. It is also recommended for infectious isolation, oncology, maternity, and Home Care. It is suitable for adults and children.
- It's unique super-absorbent pad included in each bag is designed to turn liquid wastes into a gel within seconds. Absorbs up to 300 ml of liquid wastes.

To order:

Product code: PRBV1

Case: 8 displays per case (400 bags per case)

www.hygie.com

Toll Free : 1 (866) 588-2221

Fax : (450) 444 6222

Canada, USA

GO VOM® - Vomit Bag

An effective solution to biological waste management and containment

BENEFITS


- Easy to use
- Super-Absorbent pad included in each bag that quickly absorbs up to 300 ml.
- Turn liquid wastes into a gel within seconds
- Makes the collection, transportation and disposal safer
- Limits the spread of harmful bacteria
- Blocks odors
- Reduces cleaning tasks
- Reduces risks of spilling
- Graduated lines on a transparent bag for easy measurement

IDEAL FOR

- Hospitals
- Clinics
- Home Care
- Ambulance
- Boat
- Plane
- Car

EASY TO USE

- 1) Unfold the GO VOM®
- 2) Press the edges to form a cone
- 3) Bring it close to your mouth
- 4) Close it by folding the edges of the cone and dispose.



To order:

Product code: PRBV1

Case: 8 displays per case (400 bags per case)

www.hygie.com

Toll Free : 1 (866) 588-2221

Fax : (450) 444 6222

Canada, USA